

Additional information for medical EMC of PS

The ME EQUIPMENT or ME SYSTEM is suitable for hospitals and other medical facilities (Professional healthcare / Home healthcare environment). Do not use the device near active HF SURGICAL EQUIPMENT and the RF shielded room of an ME SYSTEM, for example at magnetic resonance imaging, where the intensity of EM DISTURBANCES is high!

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally!

WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation!

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the [ME EQUIPMENT or ME SYSTEM], including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result!

© Copyright 2016 by Egston System Electronics Eggenburg GmbH
Version 01, 26.09.2016